

Product Datasheet

Goat Testosterone (T) ELISA kit (orb1671104)

Description	Goat Testosterone (T) ELISA kit
Reactivity	Goat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.1 ng/mL-20 ng/mL. Sensitivity: 0.05 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma.
Expiration Date	Please enquire.