

Product Datasheet

Goat Insulin(INS)ELISA Kit (orb1671107)

Description	Goat Insulin(INS)ELISA Kit
Reactivity	Goat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 8 μ U/mL-140 μ U/mL. Sensitivity: 5 μ U/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.