

## Product Datasheet

# Dog Transforming growth factor $\beta$ 1,TGF- $\beta$ 1 ELISA Kit (orb1671125)

|                            |                                                                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Dog Transforming growth factor $\beta$ 1,TGF- $\beta$ 1 ELISA Kit                                                                             |
| <b>Reactivity</b>          | Canine                                                                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                                                         |
| <b>Note</b>                | For research use only                                                                                                                         |
| <b>Application notes</b>   | Detect Range: 0.781 ng/mL-50 ng/mL. Sensitivity: 0.747 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 $\mu$ L. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                               |