

Product Datasheet

Chicken Tri-iodothyronine,T3 ELISA Kit (orb1671140)

Description Chicken Tri-iodothyronine,T3 ELISA Kit

Reactivity Gallus

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.5 ng/mL-8 ng/mL.
Sensitivity: 0.25 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm.
Sample Volume: 50-100 µL. Sample Type: serum, plasma, tissue homogenates.

Expiration Date Please enquire.