

Product Datasheet

Canine Tri-iodothyronine,T3 ELISA Kit (orb1671146)

Description	Canine Tri-iodothyronine,T3 ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.5 ng/mL-8 ng/mL. Sensitivity: 0.25 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma.
Expiration Date	Please enquire.