



Product Datasheet

Ferritin antibody (orb167115)


Description	Rabbit polyclonal antibody to Ferritin
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human Ferritin Light Chain (NP_000137.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	20KDa
Uniprot ID	P02792
Entrez	2512
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical analysis of formalin...



Immunofluorescence analysis of A549 cell...