

Product Datasheet

Canine Free Tri-iodothyronine Index(Free-T3) ELISA Kit (orb1671150)

Description	Canine Free Tri-iodothyronine Index(Free-T3) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 2 pmol/L-32 pmol/L. Sensitivity: 0.38 pmol/l. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma.
Expiration Date	Please enquire.