

Product Datasheet

Canine Folic acid,FA ELISA Kit (orb1671151)

Description	Canine Folic acid,FA ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.195 ng/mL-50 ng/mL. Sensitivity: 0.195 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma.
Expiration Date	Please enquire.