



Product Datasheet

HIV-1 p24 Recombinant Protein (orb1671194)


Description	HIV-1 p24 Recombinant Protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	gag
Preservatives	PBS (pH 7.4).
Form/Appearance	Liquid
Storage	Store in working aliquots at -70°C. Avoid freeze/thaw cycles. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Tag	His tag at C-terminus
Note	For research use only
Application notes	This recombinant protein can be used for WB and ELISA. For research use only.
Protein Sequence	Full length of HIV-1 Gag p24 protein
Purity	> 95% as determined by Bis-Tris PAGE
Source	E. coli
MW	24 kD
Uniprot ID	P04591
NCBI	AAB50258
Expiration Date	6 months from date of receipt.



ELISA
Validation of
HIV-1 p24
Recombinan...



SDS-PAGE
Validation of
HIV-1 p24
Recombi...



WB Validation
of HIV-1 p24
Protein with

...