


Product Datasheet

HIV-1 Envelope gp120/gp41 Recombinant Protein (orb1671263)

Description	HIV-1 Envelope gp120/gp41 Recombinant Protein
Tested Applications	CLIA, ELISA, LF, WB
Preservatives	PBS, pH 8.0, and 8 M Urea containing 0.1% ProClin300.
Form/Appearance	liquid
Concentration	batch dependent
Storage	Upon receipt, store at 2-8°C
Note	For research use only
Application notes	Suitable for use in ELISA, CLIA, Lateral Flow, Western Blot, and other immune tests
Purity	> 95% pure by SDS-PAGE. Ligand affinity purified.
Source	E. coli
MW	This recombinant protein contains immunodominant regions of HIV-1 gp120 and gp41 proteins with a MW of 47 kDa. It is recognized by positive samples in ELISA tests.
Expiration Date	6 months from date of receipt.



SDS-Page
Image