


## Product Datasheet

**SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody [T3P3F4] (orb1671387)**

## Description

SARS-CoV/SARS-CoV-2 Variant (Delta,



## Species/Host

Antibody [T3P3F4]

Camelus

## Reactivity

Virus

## Conjugation

Unconjugated

## Tested Applications

ELISA

ELISA  
Validation  
with  
Trimers  
of SARS1  
a...

## Immunogen

SARS-CoV-2 variant Spike proteins.

## Target

S

## Preservatives

SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody is supplied in PBS.

## Form/Appearance

Liquid

## Concentration

batch dependent

## Storage

SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody should be stored in working aliquots at -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

## Note

For research use only

## Isotype

sdAb

## Clonality

Recombinant

## Uniprot ID

[P0DTC2](#)

## NCBI

[QHD43419](#)

## Expiration Date

12 months from date of receipt.