

Product Datasheet

SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody [T3P1D2] (orb1671389)

Description

SARS-CoV/SARS-CoV-2 Variant (Delta,

Species/Host

Camelus

Reactivity

Virus

Conjugation

Unconjugated

Tested

Applications

ELISA

Immunogen

SARS-CoV-2 variant Spike proteins.

Target

S

Preservatives

SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer
Antibody is supplied in PBS.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer
Antibody should be stored in working aliquots at -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

sdAb

Clonality

Recombinant

Uniprot ID


[P0DTC2](#)

NCBI

[QHD43419](#)

Expiration Date

12 months from date of receipt.



ELISA
Validation
with
Trimers
of SARS1
a...