


Product Datasheet

SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody [T3P2D11] (orb1671390)

Description

SARS-CoV/SARS-CoV-2 Variant (Delta,



ELISA
Validation
with
Trimers
of SARS1
a...

	Antibody [T3P2D11]
Species/Host	Camelus
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	SARS-CoV-2 variant Spike proteins.
Target	S
Preservatives	SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody is supplied in PBS.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	SARS-CoV/SARS-CoV-2 Variant (Delta, Omicron (BA.1/BA.5)) Spike Trimer Antibody should be stored in working aliquots at -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	sdAb
Clonality	Recombinant
Uniprot ID	P0DTC2
NCBI	QHD43419
Expiration Date	12 months from date of receipt.