

Product Datasheet

**SARS-CoV/SARS-CoV-2 Variant (Beta, Delta,
Omicron (BA.5)) Spike Trimer Antibody [P1A7]
(orb1671391)**

Description

SARS-CoV/SARS-CoV-2 Variant

Spike Trimer Antibody [P1A7]

Species/Host

Camelus

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

SARS-CoV-2 variant Spike proteins.

Target

S

Preservatives

SARS-CoV/SARS-CoV-2 Variant (Beta, Delta, Omicron (BA.5))
Spike Trimer Antibody is supplied in PBS.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

SARS-CoV/SARS-CoV-2 Variant (Beta, Delta, Omicron (BA.5))
Spike Trimer Antibody should be stored in working aliquots at -20°C, stable for up to one year.
As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

sdAb

Clonality

Recombinant

Uniprot ID

P0DTC2

NCBI

QHD43419

Expiration Date

12 months from date of receipt.