



Product Datasheet

KMT5C Antibody (orb1671523)

Description	KMT5C Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing SUV420H2 protein under non-denaturing conditions.
Target	KMT5C
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q86Y97
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HEK...



Flow cytometry using the anti-SUV420H2 a...