



Product Datasheet

TDRD3 Antibody (orb1671525)


Description	TDRD3 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing TDRD3 protein under non-denaturing conditions.
Target	TDRD3
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9H7E2
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-TDRD3 anti...



Western Blot using anti-TDRD3 antibody R...



Immunofluorescence staining of HeLa cell...