



Product Datasheet

SETD7 Antibody (orb1671529)


Description	SETD7 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, FC, IF
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing SETD7 protein under non-denaturing conditions.
Target	SETD7
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q8WTS6
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-SETD7 anti...



Western Blot using anti-SETD7 antibody R...



Immunofluorescence staining of HeLa cell...