



Product Datasheet

PRDM4 Antibody (orb1671535)

Description	PRDM4 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing PRDM4 protein under non-denaturing conditions.
Target	PRDM4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9UKN5
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...



Immunofluorescence staining of HEK293 ce...