

## Product Datasheet

# H3K27me3 Antibody [BT164], Mouse IgG2b (orb1671570)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | H3K27me3 Antibody [BT164], Mouse IgG2b   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, ICC, IF, IHC, IP, WB  |
| <b>Immunogen</b>           | BT164 is a mice-derived lupus antibody that was established by selection on apoptotic chromatin. |
| <b>Preservatives</b>       | PBS with 0.02% Proclin 300.  |
| <b>Concentration</b>       | batch dependent  |
| <b>Storage</b>             | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.                 |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG2b kappa  |
| <b>Clonality</b>           | Recombinant  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |