



Product Datasheet

H3K27me3 Antibody [BT164], Rabbit IgG (orb1671571)


Description	H3K27me3 Antibody [BT164], Rabbit IgG
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC, IP, WB
Immunogen	BT164 is a mice-derived lupus antibody that was established by selection on apoptotic chromatin.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HeL...



Western Blot using anti-H3K27me3 antibod...



Flow cytometry using the anti-H3K27me3 a...