



Product Datasheet

H4 K8Ac; K12Ac; K16Ac Antibody [KM-2], Rabbit IgG (orb1671572)

Description	H4 K8Ac; K12Ac; K16Ac Antibody [KM-2], Rabbit IgG
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	KM-2 antibody was derived from an MRL/lpr lupus mouse.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Flow cytometry using the anti-H4 K8Ac; K...



Immunofluorescence staining of fixed HeL...