


Product Datasheet


CXCR1 Antibody (orb1671592)

Description	CXCR1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	This antibody was raised by immunizing BALB/c mice with a BSA-conjugated peptide comprising of the N-terminal part of IL-8RA (amino acids 1-30).
Target	CXCR1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P25024
Expiration Date	12 months from date of receipt.




Western Blot using anti-IL-8RA antibody

...



Immunofluorescence staining of fixed hum...



Flow cytometry using the Anti-IL-8RA ant...