


Product Datasheet

CCR5 Antibody (orb1671599)


Description	CCR5 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	This antibody was raised by immunizing mice with a synthetic phosphopeptide (CEAPER(pS)(pS)VYTR(pS)TGEQI(pS)VGL) corresponding to the C-terminal region of human CCR5.
Target	CCR5
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P51681
Expiration Date	12 months from date of receipt.



Western Blot using anti-CCR5 (phosphoserine...)



Immunofluorescence staining of fixed U937 cells.



Flow cytometry using the anti-CCR5 (phosphoserine...)