



Product Datasheet

HSPB2 antibody (orb167160)


Description	Rabbit polyclonal antibody to HSPB2.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human HSPB2 (NP_001532.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	20kDa
Uniprot ID	Q16082
Entrez	3316
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines.



Immunohistochemical staining of mouse heart tissue.



Immunofluorescence analysis of U2OS cells.