

Product Datasheet

TriMethyl-Histone H3-K9 antibody (orb167161)

Description Rabbit polyclonal antibody to TriMethyl-Histone H3-K9

Reactivity All, Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ChIP, DOT, ICC, IF, IHC, WB

Immunogen A synthetic trimethylated peptide around K9 of human histone H3(NP_003520.1).

Preservatives PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only


Isotype IgG

MW 17KDa


Uniprot ID **Q16695**

Dilution Range WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:50-200, ChIP: 1:50-200, ChIP


Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of 293T cell...



Dot-blot analysis of all sorts of methyl...



Chromatin immunoprecipitation analysis e...