

Product Datasheet

DiMethyl-Histone H3-K36 antibody (orb167162)

Description Rabbit polyclonal antibody to DiMethyl-Histone H3-K36

Reactivity All, Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ChIP, ICC, IF, IHC, IP, WB

Immunogen A synthetic dimethylated peptide around K36 of human histone H3(NP_003520.1).

Preservatives PBS with 0.05% proclin300,50% glycerol,pH7.3.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only


Isotype IgG

MW 17KDa


Uniprot ID **Q16695**

Dilution Range WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:50-200, ChIP: 1:20-100, ChIP


Expiration Date 12 months from date of receipt.



IHC-P staining of human lung cancer tiss...



Immunohistochemical analysis of formalin...



Immunofluorescence analysis of 293T cell...