



Product Datasheet

CD63 Antibody (orb1671621)

Description	CD63 Antibody
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	This antibody was raised by immunizing BALB/c mice with a plasma membrane preparation of the human melanoma MeWo cells.
Target	CD63
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P08962
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed A43...



Western Blot using anti-CD63 antibody NK...