

Product Datasheet

CD59 Antibody (orb1671624)

Description	CD59 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC
Immunogen	This antibody was raised by immunizing mice with stimulated human leukocytes.
Target	CD59
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgM kappa
Clonality	Recombinant
Uniprot ID	P13987
Expiration Date	12 months from date of receipt.