

Product Datasheet

CD53 Antibody (orb1671626)

Description	CD53 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IF
Immunogen	This antibody was raised by immunizing mice with stimulated human leukocytes.
Target	CD53
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2a kappa
Clonality	Recombinant
Uniprot ID	P19397
Expiration Date	12 months from date of receipt.