



Product Datasheet

POLR2L antibody (orb167163)

Description	Rabbit polyclonal antibody to POLR2L.
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human POLR2L(NP_066951.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	7kDa
Uniprot ID	P62875
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain tissue.



Immunohistochemical staining of mouse brain tissue.