


Product Datasheet

CD15 Antibody [BRA-4F1], Rabbit IgG (orb1671635)

Description	CD15 Antibody [BRA-4F1], Rabbit IgG
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP
Immunogen	This antibody was raised by immunizing mice with human CD15.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed K56...