


Product Datasheet

ITGA2B Antibody (orb1671639)

Description	ITGA2B Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	This antibody was raised by immunising mice with human CD41.
Target	ITGA2B
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P08514
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed hum...