



Product Datasheet

CD53 Antibody (orb1671659)


Description	CD53 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	This antibody was raised by immunising mice with human leukemia cells.
Target	CD53
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P19397
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-CD53 antib...



Western Blot using anti-CD53 antibody HD...



Immunofluorescence staining of fixed Dau...