

Product Datasheet

L1CAM Antibody (orb1671661)

Description	L1CAM Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	This antibody was raised by immunising mice with recombinant L1-Fc fusion protein, consisting of the ectodomain of human L1.
Target	L1CAM
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P32004
Expiration Date	12 months from date of receipt.