


Product Datasheet

Tnfsf12 Antibody (orb1671680)

Description	Tnfsf12 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	This antibody was raised by immunising SD rats with CD8-mTWEAK, a fusion of the human CD8alpha extracellular region to murine TWEAK.
Target	Tnfsf12
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	O54907
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of mouse spl...