


Product Datasheet

MS4A4A Antibody (orb1671687)

Description	MS4A4A Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	This antibody was generated by immunizing C57BL/6 mice against mouse cell line 300-19 expressing transfected human membrane-spanning 4-domains subfamily A member 4A (MS4A4A)-GFP.
Target	MS4A4A
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	Q96JQ5
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed K562...