


Product Datasheet

IL2rg Antibody (orb1671692)

Description	IL2rg Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IP
Immunogen	This antibody was raised by immunising Lewis rats with transfected rat myeloma cell line YB2/0, expressing murine cytoplasmic-tailless IL-2R common gamma-chain (YB2/0-myt).
Target	IL2rg
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P34902
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed mou...