


Product Datasheet

Cd96 Antibody (orb1671698)

Description	Cd96 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	This antibody was raised by immunising rats with murine CD96.
Target	Cd96
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	Q3U0X8
Expiration Date	12 months from date of receipt.



ELISA of anti-CD96 antibody on CD96-Fc f...