


Product Datasheet

IL17f Antibody (orb1671705)

Description	IL17f Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	This antibody was raised by immunising mice with murine IL-17F coupled to ovalbumin to generate a highly immunogenic complex.
Target	IL17f
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q7TNI7
Expiration Date	12 months from date of receipt.



Flow-cytometry using the Anti-IL-17F ant...