



Product Datasheet

Il6ra Antibody (orb1671716)

Description	Il6ra Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	This antibody was raised by immunising Lewis rats with OKT4 cells.
Target	Il6ra
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P22272
Expiration Date	12 months from date of receipt.



Flow-cytometry using the Anti-IL-6R anti...



Immunofluorescence staining of mouse spl...