


Product Datasheet

CD4 Antibody (orb1671737)

Description	CD4 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, NeA
Immunogen	This antibody was raised by immunising BALB/c mice with human CD4 complexed with HIV-1 Env (gp120).
Target	CD4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed hum...