


Product Datasheet

Cd40 Antibody (orb1671748)

Description	Cd40 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP
Immunogen	This antibody was raised by immunising Lou rats with murine CD40-hy1-Fc according to the protocol by Hasbold et al. (Eur J Immunol. 1994)
Target	Cd40
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P27512
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed RAW...