



Product Datasheet

AXIN2 antibody (orb167175)


Description	Rabbit polyclonal antibody to AXIN2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human AXIN2(NP_004646.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100KDa
Uniprot ID	Q9Y2T1
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



IHC-P staining of rat kidney tissue usin...



Immunohistochemical analysis of formalin...