

Product Datasheet

Tnfrsf9 Antibody (orb1671775)

Description	Tnfrsf9 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Activation, FC
Immunogen	This antibody was generated by immunizing rats with recombinant CD137-Fc fusion protein.
Target	Tnfrsf9
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2a kappa
Clonality	Recombinant
Uniprot ID	P20334
Expiration Date	12 months from date of receipt.