



Product Datasheet

Ly6g Antibody (orb1671788)

Description	Ly6g Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Depletion, FC, IHC, WB
Immunogen	This antibody was raised in Lou-M rats using a popliteal lymph node immunisation and fusion protocol as described by Martiniello et al (1994)
Target	Ly6g
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P35461
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of mouse spl...



Flow-cytometry using the Anti-Ly6G/Ly6C ...