



Product Datasheet

Vsir Antibody (orb1671801)

Description	Vsir Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, IHC
Immunogen	The antibody was generated by immunising Armenian hamsters with mPD-1HIg.
Target	Vsir
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	Q9D659
Expiration Date	12 months from date of receipt.



ELISA of anti-VISTA antibody on VISTA-Fc...



Immunofluorescence staining of fixed mou...