



Product Datasheet

CLEC4A Antibody (orb1671816)

Description	CLEC4A Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	This antibody was raised by immunising BALB/c mice with a dendritic cell immunoreceptor ectodomain-human Fc fusion protein.
Target	CLEC4A
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9UMR7
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of human per...



Flow-cytometry using the anti-DCIR antib...