



Product Datasheet

DARC Antibody [2C3], Rabbit IgG (orb1671828)

Description	DARC Antibody [2C3], Rabbit IgG
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC
Immunogen	CHO cells transfected with cDNA encoding the Duffy glycoprotein.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q4VBN9
Expiration Date	12 months from date of receipt.



Western Blot using anti-DARC antibody 2C...



Immunofluorescence staining of fixed Hep...