


## Product Datasheet

### Pdcd1lg2 Antibody (orb1671833)

|                            |                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Pdcd1lg2 Antibody                                                                                                           |
| <b>Reactivity</b>          | Mouse                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                |
| <b>Tested Applications</b> | Blocking, FC, IP, WB                                                                                                        |
| <b>Immunogen</b>           | This antibody was raised by immunising SD rats with B7-DC-transfected RAW264.7 cells and selected on B7-DC-transfected NRK. |
| <b>Target</b>              | Pdcd1lg2                                                                                                                    |
| <b>Preservatives</b>       | PBS with 0.02% Proclin 300.                                                                                                 |
| <b>Concentration</b>       | batch dependent                                                                                                             |
| <b>Storage</b>             | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.                                            |
| <b>Note</b>                | For research use only                                                                                                       |
| <b>Isotype</b>             | IgG kappa                                                                                                                   |
| <b>Clonality</b>           | Recombinant                                                                                                                 |
| <b>Uniprot ID</b>          | <a href="#">Q9WUL5</a>                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                             |



Binding curve of anti-PDL2 antibody TY25...