



Product Datasheet

CD74 Antibody [In-1], Rabbit IgG (orb1671845)

Description	CD74 Antibody [In-1], Rabbit IgG
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IP
Immunogen	Mouse lymphocytes of strain A injected into Wistar rats.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	CD74
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed mou...



Western Blot using anti-CD74 antibody In...